



## Nomograph


Speed of rotation (rpm)



G-force (g)



Radius of rotation



## Calculate the correct RPM

to achieve the correct g force

### Using the RCF Nomograph

To Calculate the correct RPM to achieve the correct g force, place a straightedge on the nomograph connecting the known Radius of rotation (cm) and the desired G-force (g). The point at which the straightedge intersects the Speed of rotation (rpm) axis is the rpm.

For example, if the rotating radius is 14 cm and the G-force (g) is, **300g** (Nidacon protocol), the relative centrifugal force is 1383 rpm.

### Using the formula

$$Rpm = \sqrt{\left[ \frac{g}{(1.118 \times r)} \right]} \times 10^3$$

g = the centrifugal force

r = rotational radius, the distance from the centre of the rotor to the bottom of a centrifuge tube in the bucket when raised to horizontal position

G Force /RPM calculator

<http://drycake.com/calculator/gforce.php>



[www.nidacon.com](http://www.nidacon.com)

Nidacon International AB

E-mail [contact@nidacon.com](mailto:contact@nidacon.com)

Tel: +46 31 703 06 30 Fax: +46 31 40 54 15