
PureSperm Wash™

EN

Intended Use

For washing the sperm pellet recovered from a PureSperm density gradient, for use in the swim-up method and for extending semen, or sperm pellet for use in IUI.

Components

Sodium chloride	Purified water
Potassium chloride	Glucose
Magnesium sulphate	Calcium lactate
Potassium hydrogen phosphate	Sodium pyruvate
Sodium bicarbonate	EDTA
Human serum albumin (HSA)	HEPES

Performance Characteristics

pH	7.3-8.5
Osmolality (mOsm/kg H ₂ O)	290-300
Endotoxin levels	<1.0 EU/ml
Sperm survival 18 hours after density gradient separation	>70%

Bottles and stoppers are M.E.A. tested

Storage and Stability

Store at 2-8 °C and avoid temperatures above or below these values. Under these conditions PureSperm Wash has a shelf-life of 12 months. The expiry date is shown on both bottles and cartons.

Open and close bottles under aseptic conditions. After opening, store at 2 to 8 °C when not in use. Shelf-life on the product label applies when the product is stored according to manufacturer's recommendations.

No antibiotics, additive preservatives or preservatives have been added by the manufacturer to PureSperm Wash.

Precautions and Warnings

- Use aseptic technique throughout the entire process.

- If available, use sealed buckets during centrifugation to avoid creation of aerosols.

- PureSperm Wash does not represent any fire or combustion hazard. A safety data sheet is available from the distributor or manufacturer (www.nidacom.com).

- Do not use any PureSperm Wash which shows evidence of bacterial contamination or overgrowth.

- Genital fluid may contain pathogenic bacteria. Do not touch eyes, nose or mouth.

- Do not touch eyes, nose or mouth.

- Do not re-use PureSperm Wash from any procedure due to risk for cross contamination.

- Federal Law (USA) restricts this device to sale or by or on behalf of a physician.

- Please check for regulatory requirements governing the use of ART products in your country.

- No known indication of toxic, allergic sensitivity or any pathological effects

- Standard measures to prevent infections resulting from the use of medical products prepared from human blood or plasma, including selection of donors, testing of individual donor and plasma units for specific markers of disease, and the inclusion of effective manufacturing practices for the inactivation or灭活 of viruses.

- The findes inden rapporter om behovet visivitrofertilisering med aktiverat frøststillelse i henhold til europeiske farskopspefikcerne efter fastlagte proceser.

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- There are no reports of transmission of infections with albumin manufactured to European pharmacopeia by established processes.

- It is strongly recommended that every time PureSperm Wash is administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

Gradient procedure

Recommendations

Prepare two PureSperm gradients for each semen sample. This reduces the risk of overloading a single gradient, provides security when handling tubes or recovering sperm pellets and allows two tubes to balance the centrifuge rotors.

Reagents and Equipment

- PureSperm Wash

- Dispenser sterile round bottomed centrifuge tubes

- Disposable sterile conical centrifuge tubes

- Sterile pipettes

- CO₂ incubator

- Bench top centrifuge with swing out rotor

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- Density gradient procedure

Bring all solutions to room temperature (17 - 27 °C)

1. Use a pipette with a sterile tip to add 2 mL of PureSperm 80 to a conical centrifuge tube

2. Use a new sterile pipette tip to carefully layer 2 mL PureSperm 40 on top of the PureSperm 80

3. Carefully layer liquefied semen (up to 1.5 mL) onto the PureSperm gradient

4. Centrifuge at 300 x g for 20 minutes. Do not use the brake.

5. Use a sterile Pasteur pipette to aspirate the supernatant (or the lowest 0.25 mL liquid).

6. Transfer sperm pellet to new tube and resuspend pellet in 5mL PureSperm Wash

7. Centrifuge at 500 x g for 10 minutes. Do not use the brake.

8. Carefully remove the supernatant (or the lowest 0.25 mL liquid) as possible above the pellet. If no pellet is seen, leave the bottom 0.25 mL fluid

9. Resuspend the sperm pellet in a suitable volume of culture medium to obtain the desired sperm concentration.

The sperm sample is now ready for analysis or use

To achieve the correct g force:

www.nidacom.com/rpm

Swim-up

Recommendations

Since PureSperm Wash does not contain any antibiotics, it is recommended to add antibiotics when used in preparations for ART (e.g., Penicillin, 100 U/ml).

Equalize PureSperm Wash at room temperature prior to use

Reagents and Equipment

- PureSperm Wash

- Dispenser sterile round bottomed centrifuge tubes

- Disposable sterile conical centrifuge tubes

- Sterile pipettes

- CO₂ incubator

- Bench top centrifuge with swing out rotor

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- Density gradient procedure

Bring all solutions to room temperature (17 - 27 °C)

1. Use a pipette with a sterile tip to add 2 mL of PureSperm 80 or PureSperm 100 with PureSperm Buffer

- Bench top centrifuge with swing out rotor

- Sterile pipettes

- CO₂ incubator

- Bench top centrifuge with swing out rotor

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- Density gradient procedure

Bring all the solutions to room temperature (17 - 27 °C)

1. Use a new sterile Pasteur pipette to aspirate the semen (up to 1.5 mL) onto the PureSperm gradient

2. Use a new sterile pipette tip to carefully layer 1.5 mL PureSperm Wash over the semen

3. Without disturbing the layers, place the centrifuge tube and contents, at a 45° angle, into a 5-6% CO₂ incubator at 37°C for 60 minutes.

4. Remove the lid and gently rock the centrifuge tube.

5. After centrifugation, remove the top 0.25 mL.

6. Use a new sterile Pasteur pipette to transfer the sperm pellet (or the lowest 0.25 mL liquid).

7. Resuspend the sperm pellet in a suitable volume of culture medium to obtain the desired sperm concentration.

The sperm sample is now ready for analysis or use

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