



Peroxidized mineral oil increases the oxidant status of culture media and inhibits *in vitro* porcine embryo development



C.A. Martinez ^{a, b}, A. Nohalez ^{a, b}, J.J. Ceron ^{a, b}, C.P. Rubio ^a, J. Roca ^{a, b}, C. Cuello ^{a, b},
H. Rodriguez-Martinez ^c, E.A. Martinez ^{a, b, *}, M.A. Gil ^{a, b}

^a Faculty of Veterinary Medicine, International Excellence Campus for Higher Education and Research “Campus Mare Nostrum”, University of Murcia, Murcia, Spain

^b Institute for Biomedical Research of Murcia (IMIB-Arrixaca), Murcia, Spain

^c Department of Clinical & Experimental Medicine (IKE), Linköping University, Linköping, Sweden

ARTICLE INFO

Article history:

Received 21 June 2017

Received in revised form

22 July 2017

Accepted 22 July 2017

Available online 27 July 2017

Keywords:

In vitro production of embryos

Porcine

Mineral oil overlay

Peroxidation

In vitro fertilization

In vitro maturation

ABSTRACT

The use of oils with undetected alterations is a long-recognized problem for *in vitro* embryo production systems. Since peroxides in oils have been associated with reduced embryo production outcomes, our goals were (1) to evaluate the effects of a batch of mineral oil (MO) that was suspected to be altered on the *in vitro* production of pig embryos and (2) to determine oil peroxide values throughout culture and the transfer of oxidant agents from oil to culture media. Sunflower oil, which has a completely different chemical composition than MO but a higher oxidative status, and unaltered MO were used as controls.

Oocyte maturation, fertilization and embryo development were affected differently depending on the oil overlay used. While the suspected MO was not able to sustain *in vitro* maturation and fertilization, the oocytes incubated in the presence of sunflower oil were matured and fertilized similarly to those of the unaltered MO group. Moreover, the cleavage rate of presumed zygotes cultured under the suspected MO was severely reduced compared with those cultured under the other oils, and none of the cleaved embryos developed to the blastocyst stage. Although the cleavage rates in the sunflower oil and unaltered MO groups were similar, embryos cultured under sunflower oil also failed to develop to the blastocyst stage. Our results revealed that the suspected MO and sunflower oil had similar levels of peroxides and that these levels were much higher than those of the unaltered MO. The total oxidant status was higher in media incubated under peroxidized oils than in fresh media or media incubated without an oil overlay or under unaltered MO, indicating that oxidant agents were transferred to the incubation media. However, unlike the sunflower oil group, the culture media incubated under the suspected MO had high levels of total oxidant status and low levels of hydrogen peroxide and reactive oxygen species, suggesting the presence of other unknown oxidant agents in that oil. These results indicate that a peroxidized MO overlay dramatically decreases embryo production outcomes. This decrease could be associated with the higher peroxide values of the oil but cannot be explained by the levels of hydrogen peroxide and reactive oxygen species transferred from the oil to the culture media. It is likely that different oxidant agent(s) and/or other toxic compounds present in the peroxidized MO are responsible for its damaging effects on oocytes and embryos.

© 2017 Elsevier Inc. All rights reserved.

1. Introduction

Mineral oil is widely used in embryo *in vitro* production (IVP) systems to cover drops of culture media because among other

advantages, it helps to maintain a stable and proper pH and osmolality in the culture medium [1–5]. This fact is critical since variations in the pH of the medium might affect embryo metabolism [6,7], whereas osmotic stress may cause DNA and protein alterations that affect cell function [8]. However, the use of an oil overlay may also have adverse effects due to the possibility that toxic contaminants that accumulate in the oil during production, transport or storage may be introduced into the culture medium [9,10] and may negatively affect gametes and embryos. The use of

* Corresponding author. Faculty of Veterinary Medicine, International Excellence Campus for Higher Education and Research “Campus Mare Nostrum”, University of Murcia, Murcia, Spain.

E-mail address: emilio@um.es (E.A. Martinez).

oils with undetected contamination has long been recognized as a problem for embryo cultures. Oil degradation produces alkenals and aldehydes [11], which have been shown to be toxic to embryos when present in culture media [12]. Triton X-100, a nonionic surfactant that is used in research to permeabilize membranes or solubilize proteins, has also been found in oils used for embryo culture where the embryo development was severely inhibited [11]. In another study, zinc was identified as a possible toxic contaminant of silicone oil in mouse embryo cultures [13]. However, peroxides, which produce reactive oxygen species, are possibly the most dangerous contaminants found in altered oils. They can be formed throughout the shelf life of the oil and have been associated with decreased fertilization, cleavage, and blastocyst formation rates in rodents [11,14,15]. Although these findings suggest that peroxidized oil overlays have embryo-toxic effects, the transfer of oxidant compounds from the oil to the culture medium has not been clearly established.

The aims of this study were (1) to evaluate the effects of a batch of mineral oil, which was suspected to be altered (SMO), on the IVP of pig embryos and (2) to determine oil peroxide values (POVs) throughout culture and the transfer of oxidant agents from oil to culture media.

2. Materials and methods

All experiments were performed following the ethical guidelines of the 2010/63/EU EEC Directive for animal experiments and were approved by the Ethical Committee for Experimentation with Animals of the University of Murcia, Spain (research code: 183/2015).

2.1. Collection of cumulus-oocyte complexes (COCs)

Pre-pubertal gilts of 5.5–6 months of age and weighing 90–100 kg were used as COC donors. The pre-pubertal gilts were sacrificed in a local slaughterhouse, and their ovaries were transported to the laboratory at the University of Murcia at 35 °C within 1 h post-collection in 0.9% NaCl containing 70 µg/mL kanamycin. The COCs were collected with a surgical blade from medium-sized healthy follicles (3–6 mm diameter) in Tyrode's lactate medium supplemented with 10 mM HEPES and 0.1% polyvinyl alcohol (PVA) [16,17]. Oocytes surrounded by two or more compact cumulus mass complexes and with a dark and granulated cytoplasm were selected for maturation.

2.2. In vitro maturation (IVM) and assessment of maturation

The COCs were washed three times in IVM medium. This medium consisted of TCM-199 (Gibco Life Technologies S.A. Barcelona, Spain) supplemented with 0.57 mM cysteine, 0.1% PVA and 10 ng/mL epidermal growth factor. Groups of 40 COCs were matured in 500 µL of IVM medium supplemented with 10 IU/mL equine chorionic gonadotrophin (Folligon, Intervet International B.V., Boxxmeer, The Netherlands) and 10 IU/mL human chorionic gonadotrophin (Veterin Corion, Divasa Farmavic, S.A., Barcelona, Spain) for 22 h. The COCs were then incubated for an additional 22 h in the same medium but without hormones. All cultures were incubated under an oil overlay at 38.5 °C in a humidified atmosphere of 5% CO₂ in air.

To assess the maturation status, the oocytes were fixed in acetic acid:ethanol (1:3) for 72 h at room temperature and stained with 1% lacmoid in 45% acetic acid. The oocytes were considered immature when their chromatin was enclosed in a nuclear membrane (i.e., germinal vesicle stage) or when it was condensed in metaphase I. The oocytes were considered mature when their

chromosomes were organized in metaphase and a clear first polar body was visible (MII). Oocytes with a broken cytoplasmic membrane or a cytoplasm with an abnormal appearance were considered to be degenerated. The degenerated, immature and mature oocyte rates were expressed as the ratio of the number of degenerated, immature and MII stage oocytes, respectively, relative to the total number of oocytes evaluated.

2.3. In vitro fertilization (IVF) and assessment of fertilization status

At the end of the maturation period, the COCs were denuded by vortexing at 1660 rounds/min for 2 min in IVM medium containing 0.1% hyaluronidase and washed three times in IVM medium. The oocytes were then washed three times in pre-equilibrated IVF medium [18] supplemented with 2.0 mM caffeine and 0.2% bovine serum albumin (BSA) and inseminated as described previously [19]. Briefly, the oocytes were placed in groups of 40 in 50-µL drops of IVF medium covered with oil and kept in the incubator at 38.5 °C and 5% CO₂ in humidified air for 30 min until sperm were added. Semen cryopreserved in 0.5-mL straws, as described by Carvajal et al. [20], was thawed (two straws for each replicate) at 37 °C for 20 s and washed three times at 1900 × g for 3 min in Dulbecco's PBS (Gibco, Grand Island, NY) supplemented with 0.1% BSA. The sperm pellet was then re-suspended in IVF medium, and after being appropriately diluted, 50 µL of the sperm suspension was added to the drop of IVF medium that contained the oocytes. The spermatozoa:oocyte ratio was 1000:1. The oocytes and spermatozoa were co-incubated at 38.5 °C in a humidified atmosphere of 5% CO₂ in air for 5 h under oil.

To evaluate fertilization parameters, presumptive zygotes were fixed and stained as described for the maturation assessment. Presumptive zygotes were classified as penetrated when at least a sperm head and/or male pronucleus, with the corresponding sperm tail, and two polar bodies were visible in their cytoplasm. Oocytes with an abnormal cytoplasm appearance were classified as degenerated. The penetration rate was defined as the ratio of the number of penetrated oocytes relative to the total number of mature oocytes. The monospermic rate was the ratio of the number of oocytes with one female pronucleus, one male pronucleus and two polar bodies relative to the total number of penetrated oocytes. The fertilization efficiency was defined as the ratio of the number of monospermic oocytes relative to the total number of oocytes inseminated. The degeneration rate was the ratio of the number of degenerated presumptive zygotes relative to the total number of oocytes evaluated.

2.4. In vitro culture (IVC) and assessment of embryonic development parameters

Following sperm-oocyte co-incubation, groups of 30 presumptive zygotes were washed three times and cultured in 500 µL of IVC medium, which consisted of glucose-free pre-equilibrated North Carolina State University culture medium (NCSU-23; [21]) supplemented with 0.4% BSA, 0.3 mM pyruvate and 4.5 mM lactate, for 2 days. Then, the cleaved embryos were cultured in fresh IVC medium supplemented with 0.4% BSA and 5.5 mM glucose for an additional 5 days. All embryo cultures were incubated at 38.5 °C in an atmosphere of 5% CO₂ in air under oil.

Embryonic development was evaluated morphologically under a stereomicroscope. Embryos with a morphologically abnormal appearance and those that were poorly developed were classified as degenerated embryos. The cleavage rate was defined as the ratio of the number of embryos developed to the 2- to 4-cell stage at day 2 of culture relative to the total number of oocytes inseminated. The blastocyst formation rate was the ratio of the number of 2- to 4-cell

embryos that developed to blastocyst stage (an embryo with a clear blastocoel and with good or excellent morphology) at day 7 of culture. The blastocyst efficiency rate was defined as the ratio of the number blastocysts relative to the total number of oocytes inseminated.

2.5. Analysis of oils and media

2.5.1. Determination of the peroxide values in the oils

The POVs were determined as previously reported [22]. The oil samples were treated with an acetic acid-chloroform solution and a saturated potassium iodide solution and were placed into the burette of a POV meter (DL 25 Titrator, Mettler-Toledo AG, Greifensee, Switzerland). The titration was run against a standard solution of sodium thiosulfate. The POVs were expressed in $\mu\text{mol/L}$.

2.5.2. Determination of total oxidant status (TOS), hydrogen peroxide (H_2O_2) and reactive oxygen species (ROS) in the media

The TOS values were assessed as previously described [23] using an Olympus AU600 Automatic Chemistry Analyzer. Briefly, xylenol orange, NaCl and glycerol in a H_2SO_4 solution were incubated with the sample for 3 min. Subsequently, ferrous ion and *o*-dianisidine in H_2SO_4 were added to the reaction. The change in absorbance was monitored, and the results were calculated using a standard curve of H_2O_2 solutions and expressed in $\mu\text{mol/L}$.

H_2O_2 was detected as described by Rhee et al. [24]. Briefly, the samples were incubated for 15 min at 22 °C with horseradish peroxidase (HRP) and 3,5,3',5'-tetramethylbenzidine (TBM) in sodium phosphate buffer, and the reaction was then terminated by the addition of sulfuric acid. The H_2O_2 concentration of the samples was estimated based on the molar extinction coefficient of diamine using a microplate reader (Powerwave XS, Biotek Instruments, Carson City, NV). The results were expressed in $\mu\text{mol/L}$.

ROS were measured by a luminol-mediated chemiluminescence assay [25]. Luminol containing HRP was added to the medium and saline solution. The ROS values were assessed using a microplate reader (Victor 2 1420 Multilabel Counter; PerkinElmer, Finland) and expressed in photon count per second (cps). A negative or background control was prepared by adding saline solution instead of sample.

2.6. Experimental design

In the present study, we compared the effects of an SMO overlay (Cat. no. M8410, lot no. MKBW2313V; Sigma-Aldrich Co., St. Louis, MO, USA) on IVM, IVF and embryonic development. This oil had performed poorly in previous experiments. Sunflower oil (SFO), which has a completely different chemical composition than MO but a high oxidative status, and unaltered MO (UMO; Cat. no. M8410, lot no. MKBX6122V; Sigma-Aldrich Co., St. Louis, MO, USA) were used as controls.

To evaluate the maturation and fertilization parameters, immature oocytes were matured and fertilized in media covered with SMO, SFO or UMO in four replicates. At 44 h of maturation, a random subset of the oocytes ($n = 460$) was fixed and stained to determine the percentage of oocytes reaching the MII stage. Eighteen hours post-insemination, the remaining oocytes ($n = 455$) were fixed and stained to determine the degeneration, penetration and monospermic rates, as well as the fertilization efficiency. In another set of experiments, immature oocytes ($n = 316$) were matured and fertilized under UMO and then cultured under SMO for seven days to evaluate the effects of SMO on the embryonic development. Immature oocytes matured, fertilized and cultured under SFO ($n = 305$) or UMO ($n = 413$) were used as controls.

To evaluate the oxidation status in oils and media, 500- μL drops of IVM or IVC medium covered with SMO, SFO or UMO were incubated for 22 h (IVM medium) or for 2 and 5 days (IVC medium) in the absence of oocytes or embryos. At the end of each incubation period, the oil cover was transferred to 5 mL-tubes for analysis. The medium was carefully transferred to 1 mL-tubes, avoiding cross-contamination with the oil. This process was repeated three times. The assessments were performed in three replicates immediately after sample collection.

2.7. Statistical analysis

Binary variable data (degeneration, maturation, penetration, monospermy, cleavage, blastocyst formation and efficiencies) were obtained by calculating the percentage in every drop of each group and in each of the four replicates and are expressed as the mean \pm standard deviation (SD). Continuous variables are expressed as the mean \pm SD of three replicates. Statistical analysis was performed using the IBM SPSS 19 Statistics package (SPSS, Chicago, IL, USA). Variables were analyzed using the Kolmogorov–Smirnov test for normally distributed data. Differences between means of more than two groups were determined by mixed-model analysis of variance (ANOVA) followed by the Bonferroni post hoc test. Pairwise comparisons of the means were performed using Student's *t*-tests. Differences were considered significant at $P < 0.05$.

3. Results

The maturation rates of oocytes cultured under each of these oils are presented in Fig. 1. The SMO greatly increased ($P < 0.02$) the proportion of oocytes with clear signs of degeneration at the end of the maturation period ($87.2 \pm 4.6\%$) compared with the SFO ($15.2 \pm 2.4\%$) and UMO ($8.5 \pm 2.4\%$) control groups. Consequently, there was a significant reduction ($P < 0.001$) in MII formation rates in oocytes matured in the presence of SFO ($4.2 \pm 2.2\%$) compared with oocytes matured under SFO or UMO ($74.5 \pm 7.1\%$ and $82.7 \pm 7.6\%$, respectively). The SMO overlay had no significant influence on the percentage of oocytes arrested at the germinal vesicle stage at 44 h of maturation.

All surviving oocytes from the SMO group were degenerated at 18 h post-insemination, so it was not possible to evaluate IVF parameters. By contrast, although oocytes matured and fertilized under SFO had increased ($P < 0.03$) rates of degeneration at 18 h




Fig. 1. In vitro maturation efficiency of oocytes incubated in maturation medium covered with suspected mineral oil (SMO; $n = 150$), sunflower oil (SFO; $n = 157$) or unaltered mineral oil (UMO; $n = 153$). The different letters within each variable indicate significant differences (a,b,c: $P < 0.03$; d,e: $P < 0.001$). The data are presented as the mean \pm SD of four replicates.




Fig. 2. In vitro fertilization parameters of oocytes matured and fertilized in media covered with suspected mineral oil (SMO; n = 156), sunflower oil (SFO; n = 152) or unaltered mineral oil (UMO; n = 147). Penetration: percentage of oocytes penetrated/total inseminated oocytes. Monospermy: percentage of oocytes containing one female pronucleus, one male pronucleus and two polar bodies/total penetrated oocytes. Fertilization efficiency: percentage of monospermic oocytes/total inseminated oocytes. The different letters within each variable indicate significant differences ($P < 0.001$). The data are presented as the mean \pm SD of four replicates.

post-insemination compared with the UMO group ($13.2 \pm 5.4\%$ vs. $5.0 \pm 2.2\%$), no differences were observed in rates of sperm penetration, monospermy and IVF efficiency between these two groups (Fig. 2).

To evaluate the effects of SFO on embryonic development, immature oocytes were matured and fertilized under UMO and then cultured for 7 days under SMO. Embryonic development data are presented in Fig. 3. There was no difference in cleavage rate between the SFO and UMO groups ($62.9 \pm 7.8\%$ and $59.6 \pm 6.9\%$, respectively), but the proportion of cleaved embryos derived from presumed zygotes cultured under SMO ($25.0 \pm 4.9\%$) was dramatically decreased ($P < 0.001$). The blastocyst formation rates achieved in the UMO group ($33.1 \pm 7.3\%$) were within the normal range in our laboratory. However, no blastocyst formation was observed upon culture of presumed zygotes under SMO or SFO.

The POVs in oils incubated with IVM or IVC medium are presented in Tables 1 and 2. At 0 h of incubation, the POVs of SMO and SFO were approximately 100-fold greater than those of UMO. These differences in POVs were maintained throughout the incubation period. Although no differences in POVs were observed between




Fig. 3. Embryonic development of presumed zygotes cultured for seven days in culture medium covered with suspected mineral oil (SMO; n = 316), sunflower oil (SFO; n = 305) or unaltered mineral oil (UMO; n = 413). Presumed zygotes from the SMO group were matured and fertilized in media covered with UMO. Blastocyst formation: percentage of blastocysts/total cleaved embryos. Blastocyst efficiency: percentage of blastocysts/total cultivated oocytes. The different letters within each variable indicate significant differences ($P < 0.001$). The data are presented as the mean \pm SD of four replicates.

Table 1
Peroxide values ($\mu\text{mol/L}$) in oils incubated with maturation medium for 22 h.

Group	Time of culture (h)	
	0	22
Suspected MO	$638.4 \pm 75.6^{\text{aA}}$	$803.7 \pm 55.6^{\text{aB}}$
Sunflower oil	$565.7 \pm 87.7^{\text{a}}$	$626.1 \pm 35.2^{\text{b}}$
Unaltered MO	$6.4 \pm 0.6^{\text{b}}$	$6.1 \pm 0.1^{\text{c}}$

MO: mineral oil.

^{a,b,c}Different letters in the same column indicate differences ($P < 0.001$).

^{A,B}Different letters in the same row indicate differences ($P < 0.05$).

Data are presented as the mean \pm SD (three replicates).

Table 2
Peroxide values ($\mu\text{mol/L}$) in oils incubated with embryo culture medium for 2 and 5 days.

Group	Time of culture (days)		
	0	2	5
Suspected MO	$638.4 \pm 75.6^{\text{a}}$	$820.3 \pm 106.0^{\text{a}}$	$789.6 \pm 73.2^{\text{a}}$
Sunflower oil	$565.7 \pm 87.7^{\text{a}}$	$592.7 \pm 82.3^{\text{b}}$	$585.0 \pm 28.1^{\text{b}}$
Unaltered MO	$6.4 \pm 0.6^{\text{b}}$	$6.8 \pm 0.7^{\text{c}}$	$6.6 \pm 0.4^{\text{c}}$

MO: mineral oil.

^{a,b,c}Different letters in the same column indicate differences ($P < 0.008$).

Data are presented as the mean \pm SD (three replicates).

SMO and SFO at 0 h of incubation, the POVs were close to 1.5-fold higher ($P < 0.03$) in SMO than in SFO after incubation, regardless of culture time and medium.

The oxidant levels of IVM and IVC media incubated under the different oils are presented in Tables 3 and 4. When we evaluated the oxidative status of fresh IVM and IVC media (0 h of incubation), the levels of TOS and H_2O_2 were below the limit of detection. These values remained below the limit of detection at 22 h of incubation (IVM medium), while the TOS levels increased above this limit at day 2 or 5 of incubation (IVC medium). The ROS values in fresh IVM medium (4872.5 ± 57.6 cps) increased ($P < 0.05$) to 5046.7 ± 167.6 cps at 22 h of incubation in the absence of an oil overlay. A similar pattern was observed for IVC medium, in which the ROS values increased ($P < 0.001$) from 3447.0 ± 40.8 cps at 0 h of incubation to 5483.3 ± 90.1 cps at 2 days of incubation in the absence of an oil overlay. The TOS, H_2O_2 and ROS values of IVM and IVC media incubated under UMO were similar to those observed in the media incubated without an oil overlay, regardless of the incubation time. By contrast, IVM medium or IVC medium incubated under SMO had higher values ($P < 0.02$) of TOS than media incubated under UMO, while the H_2O_2 and ROS values were not significantly different. Media incubated under SFO had higher ($P < 0.02$) levels of TOS, H_2O_2 and ROS than media incubated without an oil overlay or medium covered with UMO at any incubation time.

4. Discussion

The results from this study clearly demonstrate that porcine oocytes and embryos did not survive IVC due to one batch of mineral oil that had been used as an oil overlay for *in vitro* maturation, fertilization and embryo culture. This study also supports previous findings in mice regarding the negative effects of peroxide-contaminated oil on *in vitro* fertilization and embryo development success [14]. In addition, we also demonstrate the transfer of oxidant molecules from the oil overlay into the incubation media.

This study was conceived because in our routine experiments, we noticed a sudden decrease in maturation, fertilization and embryonic development rates compared with our previous

Table 3

Total oxidant status (TOS; $\mu\text{mol/L}$), hydrogen peroxide (H_2O_2 ; $\mu\text{mol/L}$) and reactive oxygen species (ROS; cps) levels in maturation media covered with the different oils and incubated for 22 h.

Oil cover	Time of culture (hours)					
	0			22		
	TOS	H_2O_2	ROS	TOS	H_2O_2	ROS
Suspected MO	–	–	–	9.4 ± 0.8^a	<LOD	5191.3 ± 42.4^a
Sunflower oil	–	–	–	21.6 ± 0.6^b	9.6 ± 0.1	7264.6 ± 150.0^b
Unaltered MO	–	–	–	<LOD	<LOD	5149.4 ± 83.7^a
None	<LOD	<LOD	4872.5 ± 57.6^A	<LOD	<LOD	5046.7 ± 167.6^{AB}

MO: mineral oil; LOD: Limit of detection ($1.41 \mu\text{mol/L}$ for TOS and $1.79 \mu\text{mol/L}$ for H_2O_2).

^{a,b}Different letters in the same column indicate differences ($P < 0.001$).

^{A,B}Different letters in the same row indicate differences ($P < 0.05$).

Data are presented as the mean \pm SD (three replicates).

Table 4

Total oxidant status (TOS; $\mu\text{mol/L}$), hydrogen peroxide (H_2O_2 ; $\mu\text{mol/L}$) and reactive oxygen species (ROS; cps) levels in embryo culture media covered with the different oils and incubated for 2 and 5 days.

Oil cover	Time of culture (days)								
	0			2			5		
	TOS	H_2O_2	ROS	TOS	H_2O_2	ROS	TOS	H_2O_2	ROS
Suspected MO	–	–	–	16.6 ± 0.2^a	<LOD	5622.7 ± 231.0^a	15.7 ± 0.4^a	<LOD	5846.0 ± 574.1^a
Sunflower oil	–	–	–	4.6 ± 0.4^b	4.2 ± 0.3	16481.7 ± 743.2^b	4.3 ± 0.2^b	4.3 ± 0.4	15639.3 ± 705.0^b
Unaltered MO	–	–	–	2.8 ± 0.1^c	<LOD	5831.0 ± 156.4^a	2.9 ± 0.1^c	<LOD	6057.7 ± 194.4^a
None	<LOD	<LOD	3447.0 ± 40.8^A	2.3 ± 0.1^c	<LOD	5483.3 ± 90.1^{AB}	2.3 ± 0.3^c	<LOD	5544.7 ± 80.2^{AB}

MO: mineral oil; LOD: Limit of detection ($1.41 \mu\text{mol/L}$ for TOS and $1.79 \mu\text{mol/L}$ for H_2O_2).

^{a,b}Different letters in the same column indicate differences ($P < 0.02$).

^{A,B}Different letters in the same row indicate differences ($P < 0.001$).

Data are presented as the mean \pm SD (three replicates).

outcomes. This problem became more pronounced over time. We suspected that the problem could be related to the use of a new batch of mineral oil, since no modifications in laboratory environmental conditions or methodological procedures had been introduced. Because peroxides are one of the most damaging contaminants found in oil [11,14,15], we decided to investigate the POV of the SMO and the possible transfer of oxidant agents from the oil into the incubation media. Because gametes and embryos produce many oxidant agents via their metabolic processes (reviewed by Takahashi [26]), in our study, the incubation of medium under oil was performed in the absence of oocytes or embryos. Therefore, the oxidative status of each medium was determined in the absence of metabolic processes. Additionally, we used SFO and UMO for comparison. We used SFO for two main reasons: First, SFO, which contains mainly saturated, monounsaturated and polyunsaturated fatty acids [27,28], has a completely different chemical composition than mineral oil; and second, the POVs of SFO are usually much higher [27,29] than those of UMO.

Our results revealed that SMO and SFO had similar levels of peroxides at 0 h of incubation and that in both cases, these levels were much higher than those found in UMO. Additionally, the TOS values were higher in the media incubated under peroxidized oils compared with those obtained in fresh media (0 h of incubation) or media incubated without oil overlay or under UMO, which indicates a clear transfer of oxidant molecules from the oils into the incubation media.

Despite these findings, SMO and SFO affected oocyte maturation, fertilization and development differently. Thus, while SMO was not able to sustain *in vitro* maturation and fertilization, with all oocytes being degenerated at 18 h post-insemination, only 15% of oocytes degenerated in the presence of SFO, and the remaining oocytes were fertilized similarly to those of the UMO group. Moreover, the cleavage rates of presumed zygotes cultured under SMO were severely reduced compared with the SFO group, and

none of the cleaved embryos developed to the blastocyst stage in the SMO group. The harmful effects of SFO on embryonic development were not observed at 48 h of culture, since the cleavage rates in that group were no different from those observed in the UMO group. However, like in the SMO group, all embryos cultured under SFO failed to develop to the blastocyst stage.

The reasons for these differences between SMO and SFO are not clear. It is possible that the increase observed in the POV of SMO both at 22 h of maturation and at day 2 of culture compared with that of SFO could be responsible for the differences in degeneration and cleavage rates, respectively, between the oils. Supporting this speculation, it has been demonstrated that embryonic degeneration is accelerated with increased levels of oil peroxidation [14]. However, this is not necessarily the case, since the IVM or IVC medium incubated under SMO had similar levels of H_2O_2 and ROS as the UMO group and lower levels than the SFO group. Moreover, despite the huge differences in POV, H_2O_2 and ROS levels between the SFO and UMO groups, no differences in maturation rate, fertilization parameters or cleavage rate were observed between these groups. These findings indicate that H_2O_2 and ROS were not directly responsible for the destructive effects of SMO on oocytes and embryos.

Interestingly, the media incubated under SMO had higher levels of TOS than the media incubated under UMO. However, the results obtained regarding TOS levels in IVM and IVC media incubated under SMO seem contradictory, because these values were lower at 22 h of maturation and higher at 2 days and 5 days of embryo culture than those obtained in the SFO group. This apparent discrepancy can be explained by two reasons: First, although the main TOS contributors are hydrogen peroxide and lipid hydroperoxides [23], other components that contribute to TOS values remain unknown, which could explain the high TOS levels in the presence of low levels of H_2O_2 and ROS in the SMO group; and second, some of these unknown compounds could also interact differently

depending on the medium composition, resulting in the differences in TOS values observed in the media in this study. This hypothesis is supported by previous findings reporting that the generation of oxidant agents depends on the composition of the culture medium [30–32] and that specific substances, including ionic elements and serum albumin, present in the culture medium can increase the oxidative status of the medium [33,34].

Overall, our results suggest that the negative effects of SFO could be due, at least in part, to the increased H₂O₂ and ROS values in the culture media, which are probably associated with the higher POVs conserved in SFO. By contrast, unknown oxidant compounds detected by TOS analysis could be responsible for the poor performance of SMO. Nevertheless, the possibility that other toxic compounds released from SMO could be implicated in its damaging effects should be considered. In accordance with this hypothesis, aldehydes and alkenals, which are embryo toxic [12], have also been identified in peroxidized oils [11].

5. Conclusions

In conclusion, this manuscript reports a case in which alterations in the state of mineral oil used to overlay culture media dramatically decreased IVP outcomes. This decrease could be related to the higher POV concentration of the suspected oil but cannot be explained by the levels of H₂O₂ and ROS in the culture media. It is likely that different oxidant agent(s) and/or other toxic compounds present in the peroxidized mineral oil are responsible for its damaging effects on oocytes and embryos. Appropriate conservation and control procedures, including measuring the POV in the oils used for IVP culture media, are recommended.

Acknowledgements

The authors are grateful to AIM Iberica (Murcia, Spain) and El Pozo (Murcia, Spain) for supplying the boar ejaculates and the ovaries, respectively, used in this study.

Funding

This study was supported by the Ministry of Economy and Competitiveness (Madrid, Spain) and the European Regional Development Fund (grant number AGL2015-69735-R); and the Seneca Foundation, Murcia, Spain (grant number 19892/GERM/15). The Ministry of Economy and Competitiveness is also acknowledged for its grant-based support to C.A. Martinez and A. Nohalez (BES-2013-064087 and BES-2013-064069, respectively).

Role of the funding source

Funding sources did not have any involvement in the study design, in the collection, analysis and interpretation of data, in the writing of the report, and in the decision to submit the article for publication.

Author contributions

C.A.M., H.R.-M., E.A.M. and M.A.G. conceived the study. C.A.M., A.N., I.P., J.R., C.C., E.A.M. and M.A.G. performed the experiments related to embryo production. J.J.C. and C.P.R. collected and interpreted the data regarding oil and medium analyses. C.A.M., H.R.-M., E.A.M. and M.A.G. wrote the manuscript. All authors reviewed the manuscript.

Declaration of interest

None of the authors have any conflicts of interest to declare.

References

- Brinster RL. *In vitro* culture of mammalian embryos. *J Anim Sci* 1968;27:1–14.
- Boatman DE. *In vitro* growth of non-human primate pre- and peri-implantation embryos. In: Bavister BD, editor. The mammalian preimplantation embryo: regulation of growth and differentiation *in vitro*. New York: Plenum Press; 2009. p. 273–303.
- Miller KF, Goldberg JM, Collins RL. Covering embryo cultures with mineral oil alters embryo growth by acting as a sink for an embryotoxic substance. *J Assist Reprod Genet* 1994;11:342–5.
- Gardner DK, Lane M. *Ex vivo* early embryo development and effects on gene expression and imprinting. *Reprod Fertil Dev* 2005;17:361–70.
- Tae JC, Kim EY, Lee WD, Park SP, Lim JH. Sterile filtered paraffin oil supports *in vitro* developmental competence in bovine embryos comparable to co-culture. *J Assist Reprod Genet* 2006;23:121–7.
- Edwards LJ, Williams DA, Gardner DK. Intracellular pH of the preimplantation mouse embryo: effects of extracellular pH and weak acids. *Mol Reprod Dev* 1998;50:434–42.
- Lane M, Gardner DK. Regulation of ionic homeostasis by mammalian embryos. *Semin Reprod Med* 2000;18:195–204.
- Burg MB, Ferraris JD, Dmitrieva NI. Cellular response to hyperosmotic stresses. *Physiol Rev* 2007;87:1441–74.
- Lee S, Cho M, Kim E, Kim T, Lee C, Han J, et al. Renovation of a drop embryo cultures system by using refined mineral oil and the effect of glucose and/or hemoglobin added to a serum-free medium. *J Vet Med Sci* 2004;66:63–6.
- Morbeck DE, Leonard PH. Culture systems: mineral oil overlay. *Methods Mol Biol* 2012;912:325–31.
- Morbeck DE, Khan Z, Barnidge DR, Walker DL. Washing mineral oil reduces contaminants and embryotoxicity. *Fertil Steril* 2010;94:2747–52.
- Hall J, Gilligan A, Schimmel T, Cecchi M, Cohen J. The origin, effects and control of air pollution in laboratories used for human embryo culture. *Hum Reprod* 1998;13:146–55.
- Erbach GT, Bhatnagar P, Baltz JM, Biggers JD. Zinc is a possible toxic contaminant of silicone oil in microdrop cultures of preimplantation mouse embryos. *Hum Reprod* 1995;10:3248–54.
- Otsuki J, Nagai Y, Chiba K. Peroxidation of mineral oil used in droplet culture is detrimental to fertilization and embryo development. *Fertil Steril* 2007;88:741–3.
- Otsuki J, Nagai Y, Chiba K. Damage of embryo development caused by peroxidized mineral oil and its association with albumin in culture. *Fertil Steril* 2009;91:1745–9.
- Funahashi H, Ekwall H, Rodriguez-Martinez H. Zona reaction in porcine oocytes fertilized *in vivo* and *in vitro* as seen with scanning electron microscopy. *Biol Reprod* 2000;63:1437–42.
- Martinez EA, Angel MA, Cuello C, Sanchez-Orsorio J, Gomis J, Parrilla I, et al. Successful non-surgical deep uterine transfer of porcine morulae after 24 hour culture in a chemically defined medium. *PLoS One* 2014;9:e104696.
- Abeydeera LR, Day BN. *In vitro* penetration of pig oocytes in a modified Tris-buffered medium: effect of BSA, caffeine and calcium. *Theriogenology* 1997;48:537–44.
- Gil MA, Abeydeera LR, Day BN, Vazquez JM, Roca J, Martinez EA. Effect of the volume of medium and number of oocytes during *in vitro* fertilization on embryo development in pigs. *Theriogenology* 2003;60:767–76.
- Carvajal G, Cuello C, Ruiz M, Vázquez JM, Martínez EA, Roca J. Effects of centrifugation before freezing on boar sperm cryosurvival. *J Androl* 2004;25:389–96.
- Petters RM, Wells KD. Culture of pig embryos. *J Reprod Fertil Suppl* 1993;48:61–73.
- Association of Analytical Communities. Official methods of analysis of association of official analytical chemists. Arlington, VA, USA: Association of Analytical Communities; 1995.
- Erel O. A new automated colorimetric method for measuring total oxidant status. *Clin Biochem* 2005;38:1103–11.
- Rhee SG, Chang TS, Jeong W, Kang D. Methods for detection and measurement of hydrogen peroxide inside and outside of cells. *Mol Cells* 2010;29:539–49.
- Vong L, Lorentz RJ, Assa A, Glogauer M, Sherman PM. Probiotic *Lactobacillus rhamnosus* inhibits the formation of neutrophil extracellular traps. *J Immunol* 2014;192:1870–7.
- Takahashi M. Oxidative stress and redox regulation on *in vitro* development of mammalian embryos. *J Reprod Dev* 2012;58:1–9.
- Pal US, Patra RK, Sahoo NR, Bakhara CK, Panda MK. Effect of refining on quality and composition of sunflower oil. *J Food Sci Technol* 2015;52:4613–8.
- Orsavova J, Misurcova L, Ambrozova JV, Vicha R, Mlcek J. Fatty acids composition of vegetable oils and its contribution to dietary energy intake and dependence of cardiovascular mortality on dietary intake of fatty acids. *Int J Mol Sci* 2015;16:12871–90.
- Crapiste GH, Brevedan MIV, Carelli AA. Oxidation of sunflower oil during storage. *J Am Oil Chem Soc* 1999;76:1437.
- Agarwal A, Said TM, Bedaiwy MA, Banerjee J, Alvarez JG. Oxidative stress in an

- assisted reproductive techniques setting. *Fertil Steril* 2006;86:503–12.
- [31] Martin-Romero FJ, Miguel-Lasobras EM, Dominguez-Arroyo JA, Gonzalez-Carrera E, Alvarez IS. Contribution of culture media to oxidative stress and its effect on human oocytes. *Reprod Biomed Online* 2008;17:652–61.
- [32] Shih Y-F, Lee T-H, Liu C-H, Tsao H-M, Huang C-C, Lee M-S. Effects of reactive oxygen species levels in prepared culture media on embryo development: a comparison of two media. *Taiwan J Obstet Gynecol* 2014;53:504–8.
- [33] Shannon P. Factors affecting semen preservation and conception rates in cattle. *J Reprod Fertil* 1978;54:519–27.
- [34] Nasr-Esfahani M, Johnson MH, Aitken RJ. The effect of iron and iron chelators on the *in-vitro* block to development of the mouse preimplantation embryo: BAT6 a new medium for improved culture of mouse embryos *in vitro*. *Hum Reprod* 1990;5:997–1003.