

# QUALITY ASSURANCE CERTIFICATE



NidOil™

An oil overlay for use during ICSI or gamete and embryo culture, inside or outside the incubator.

Containing: 4x100 mL

Catalogue no: NO-400K

Batch:  
**100NOGB10**

Production date:  
**10 February 2026**

Expiry date:  
**09 February 2028**

We certify that this product was manufactured according to set directives for medical devices. The quality of this particular batch was tested and evaluated objectively by stringent, quality control procedures.

**Quality Control procedures:**

**Accepted values:**

**Test results:**

**Physical analyses**

Density (d 20°C / 20°C, g/mL)

0.82-0.88

0.849

**Sterility and toxin analyses**

Microbiological growth control

No growth

Comply

Endotoxins, quantitative LAL assay (EU/mL)

< 0.100

<0.005

**Biological analyses**

Mouse Embryo Assay

-development to blastocyst stage at day 5 (%)

>80

Comply

Human sperm survival after 18 hrs

>85

Comply

**Peroxide analyses**

Peroxide level (mEq/kg)

<0.02

0.0012

**Additional information**

Extended Mouse Embryo Assay (eMEA), development to blastocyst stage at day 6

95%

Date of release: 20 March 2026

This document is electronically generated and does not bear a signature.

**Marina Danilova**

Quality Assurance Coordinator

Nidacon International AB  
Flöjelbergsgatan 16B, SE-431 37 Mölndal, Sweden  
Tel: +46-31-703 06 30 Fax: +46-31-40 54 15  
E-mail: [contact@nidacon.com](mailto:contact@nidacon.com) Web: [www.nidacon.com](http://www.nidacon.com)